Total solution of nucleic acid chemistry by LGC, Biosearch Technologies™

Your partner of choice for all your
oligonucleotide synthesis reagent
and instrument needs.

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For Research Use Only.
Not for use in diagnostic procedures.

Catalogue

Who we are	2
Quality at Biosearch Technologies	2
Enhanced specification	3
Custom synthesis and bulk packaging	3
Chemistry of oligonucleotide synthesis	4
Choosing a solid support	4
CPG products	5
Employing universal supports	6
Bespoke 3'-incorporation: employing unfunctionalized Amino-CPG	6
Applications of modified oligonucleotides	7
Therapeutics	7
Antisense therapy	7
siRNA	8
Aptamers	8
CRISPR/Cas gene editing	8
Modifiers and their use in oligonucleotide synthesis	9
Backbone modification-Peptide Nucleic Acid (PNA)	9
Locked nucleic acid (LNA™) oligonucleotides	9
Alpha-L-LNA	10
Phosphorothioates	10
Modification for nuclease resistance	11
2'-O-Methyl modification	11
2'-Fluoro modification	12
Reverse (5' to 3') oligonucleotides	12
Cell delivery and uptake	13
Lipophilic modification	13
GalNac	14
DNA/RNA Synthesizers	20
Miscellaneous products	21

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Table 1. Summary of typical CPGs available.

Type	Modification
CPG 500 Å, 600 Å, 1000 Å	Native, AMP, CNA/LCAA
CPG 2000 Å, 3000 Å	Native, AMP, CNA/LCAA
DNA CPG 500 Å, 600 Å, 1000 Å	dA, dC (Ac), dC (Bz), dG (iBu), dG (dmf), dT
DNA CPG 2000 Å, 3000 Å	dA, dC (Ac), dC (Bz), dG (iBu), dG (dmf), dT
RNA CPG 1000 Å	rA (Bz), rC (Ac), rC (Bz), rG (iBu), rG (dmf), rU, rA (Bz), 2'-OAc, rA, (Ac), 2'-OAc, rC (Ac), 2'-OAc, rG (iBu), 2'-OAc, rU, 2'-OAc
RNA CPG 2000 Å, 3000 Å	rA (Bz), rC (Ac), rC (Bz), rG (iBu), rG (dmf), rU, rA (Bz), 2'-OAc, rA (Ac), 2'-OAc, rC (Ac), 2'-OAc, rG (iBu), 2'-OAc, rU, 2'-OAc
RNA CPG (Pac); all pore sizes	rA (Pac), rG (iPr-Pac), rG (iPr-Pac), 2'-OAc
2'-OME CPG 1000 Å	2' OMe rA (Bz), 2' OMe rC (Ac), 2' OMe rC (Bz), 2' OMe rG (iBu) 2' OMe rG (dmf), 2' OMe rU
2'-OME CPG 2000 Å, 3000 Å	2' OMe rA (Bz), 2' OMe rC (Ac), 2' OMe rC (Bz), 2' OMe rG (iBu), 2' OMe rG (dmf), 2' OMe rU
Reverse CPG	3' DMT dT, 3' DMT dA (Bz), 3' DMT dC (Bz)
LNA CPG	LNA A (Bz), LNA C (Bz), LNA G (dmf), LNA T
2'-F CPG	2' Fluoro A (Bz), 2' Fluoro C (Ac), 2' Fluoro G (iBu), 2' Fluoro U

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

MerMade™ 6 & 12

- 50 nmol to 200 μ mol+ scale synthesis
- 6 or 12 columns per run
- Reagent calibration by curve
- Automatic wobble delivery
- Up to 24 amidite positions
- In-line trityl monitoring
- Sequence manager software
- Low operational cost per oligo

Injection head


Column detail

MerMade 4

- 4 column synthesizer
- 20 nmol to 5 μ mol scale synthesis
- 10 amidite positions
- Low reagent consumption
- Fast synthesis cycle: 20mer primers in 3 hours
- In-line trityl monitoring
- Low operational cost per oligo
- Updated user friendly software

Column detail


Bottle detail

MerMade 48X

- 48 column synthesizer (4 banks of 12)
- 20 nmol to 5 μ mol scale synthesis
- Low reagent consumption
- Fast synthesis cycle: 20mer primers in 2.5hours
- User friendly control software
- Up to 24 amidite positions
- Small footprint
(D: 59 cm; W: 71 cm; H: 48 cm)
- Low operational cost per oligo
- Automatic wobble delivery

Columns

MerMade 192R

- 192 column synthesizer
- 50 nmol to 1 μ mol scale synthesis
- 6 amidite positions
- Robot loadable
- Fast synthesis cycle: 20mer primers in 3 hours
- LIMS system compatible for post processing
- Fully featured API for remote monitoring and control

MerMade 192E

Top view (injection head detail)

- 192 column synthesizer
 - 25 nmol to 1 μ mol scale synthesis (up to 5 μ mol on modified MM192Es)
 - Fast synthesis cycle: 20mer primers in 3 hours
 - Reagent calibration by curve
 - Automatic wobble delivery
 - 8 amidite positions
 - Low volume dispense option
 - Groups of 8 reagent injection
- 96E also available as a single plate with 12 amidite positions

MerMade 192X

- 192 column synthesizer
- 50 nmol to 1 μ mol scale synthesis (up to 5 μ mol on modified MM192X)
- Up to 64 amidite positions
- Argon plate pressurise feature
- Proprietary 8 channel reagent delivery system - minimises valves
- User friendly control software
- 6 ancillary + 2 auxiliary ports
- Injection head customisable - group reagent injections according to needs

Injection head detail

Amidite bottles

Table 6. Summary of automated DNA/RNA MerMade synthesizer range from Biosearch Technologies. This table is for reference purposes only and should not be used for instrument buying decisions. Please check with the Customer Service for detailed and current specification information as this may change.

Model	Scale	Format	Trityl monitor	Amidite ports (A, G, C, T/U) ports	Amidite ports (modifiers)	Ancillary reagent ports	Standard column type
MerMade 6	50 nmol - 200 µmol	Column	Yes	10 to 20		7	Pipette-tip
MerMade 12	50 nmol - 200 µmol	Column	Yes	10 to 20		7	Pipette-tip
MerMade 48X	50 nmol - 1 µmol	Column	No	2 x 4	6	7	Pipette-tip
MerMade 192X	5 nmol - 1 µmol	Plate or Column	No	10 to 20		7	Pipette-tip (if column)
MerMade 192E	50 nmol - 1 µmol	Column	No	4		7	Pipette-tip

Standard bottle type (amidites)	Standard bottle type (ancillary reagents)	Synthesis positions	Notes and availability
45 mm (GL45) or 28-405 screw	Activator/Caps/Oxidiser: 28-405 screw; Wash/Deblock: 45 mm (GL45) screw	6 expandable to 12	Dedicated sulphurisation port. For availability contact Customer Service
45 mm (GL45) or 28-405 screw	Activator/Caps/Oxidiser: 28-405 screw; Wash/Deblock: 45 mm (GL45) screw	12	Dedicated sulphurisation port. For availability contact Customer Service
20 mm slider	Activator/Caps/Oxidiser: 28-405 screw; Wash/Deblock: 45 mm (GL45) screw	48	For availability contact Customer Service
45 mm (GL45) or 28-405 screw	Activator/Caps/Oxidiser: 28-405 screw; Wash/Deblock: 45 mm (GL45) screw	2 x 96	Dedicated sulphurisation port. For availability contact Customer Service
45 mm (GL45) or 28-405 screw	Activator/Caps/Oxidiser: 28-405 screw; Wash/ Deblock: 45 mm (GL45) screw	192	For availability contact Customer Service

Miscellaneous products

We also provide nucleoside synthesis reagents, Molecular Traps, plus empty synthesizer bottles and columns.

Integrated tools. Accelerated science.

   @LGCBiosearch

bioautomation.com
biosearchtech.com

All trademarks and registered trademarks mentioned herein are the property of their respective owners. All other trademarks and registered trademarks are the property of LGC and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying, recording or any retrieval system, without the written permission of the copyright holder. © LGC Limited, 2022. All rights reserved.

BIOSEARCH™
TECHNOLOGIES
GENOMIC ANALYSIS BY LGC

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Shigematsu & Co., Ltd.