SmartEnzymes™

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Nature offers a vast source of enzymes, perfected through evolution to perform defined reactions.

At Genovis, we believe that enzymes with unique properties can be used as biological tools to support the research and development of complex biopharmaceuticals to help bring safe and effective medicines to patients in need.

Our task is to identify new enzymes and give them names.

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

Determining antibody quality attributes. Relative abundance of glycans on cetuximab Fc/2 fragments, generated with FabRICATOR and analyzed with LC-MS at the subunit level by Ayoub and colleagues (Ayoub, 2013).

Adalimumab (IgG1)


Quality control of antibody oxidation. Oxidation analysis of adalimumab showing UV (280 nm) chromatograms before (control) and after (oxidized) incubation with 0.25% H₂O₂. A deconvoluted mass spectrum of the deglycosylated Fc/2 fragments demonstrates the degree of oxidation.

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References

Ayoub, D. et al., 2013. *Correct primary structure assessment and extensive glyco-profiling of cetuximab by a combination of intact, middle-up, middle-down and bottom-up ESI and MALDI mass spectrometry techniques.* mAbs, 5(5), pp.699–710.

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


FabRICATOR-HPLC in a middle-level workflow. Potential set-up for an automated middle-level workflow using FabRICATOR-HPLC that digests IgG to F(ab')₂ and Fc, which is well suited for high resolution MS analysis.

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

Performance testing of FabRICATOR-HPLC.

Comparison between digestion with standard in-solution FabRICATOR protocol (orange) and digestion using the automated FabRICATOR-HPLC workflow (blue), and quantification of the column's performance over a period of 14 days (bottom). The asterisk marks LC fragments that are not completely reduced with one intramolecular disulfide bridge intact.

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Glycan analysis using FabRICATOR-HPLC.

Deconvoluted mass spectra of trastuzumab Fc/2 from in-solution FabRICATOR digestion (top, orange), FabRICATOR-HPLC (top, blue), and glycosylation profiles of trastuzumab (bottom, orange, n=10) or (bottom, blue, n=28).

FabRICATOR®-HPLC

Column hardware: PEEK/biocompatible

Column dimensions: 2.1 mmD × 50 mL

Support resin: POROS®*

Typical flow rate: 0.025-0.05 mL/min

Maximum Pressure: 100 bar

Operating pH: 6.5 - 8.0

Operating temperature: 37°C

Storage conditions: +4-8°C
(Do not freeze!)

Number of days of continuous operation: >10**

Injections per column: >200**

Start material: Human IgG1-4,
Fc-fusion proteins

* See legal and disclaimers. ** Depending on specific application

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

Generation of antibody subunits with FabRICATOR MagIC. Analysis of a monoclonal antibody by RPLC-MS at the intact level (top), after FabRICATOR MagIC digestion (middle) and after treatment with FabRICATOR MagIC and 5 mM TCEP (bottom). Data was obtained in collaboration with Symphogen A/S (Ballerup, Denmark).

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

Oxidation studies using FabRICATOR MagIC. Six mAbs were analyzed for tryptophan oxidation during forced degradation studies by RPLC-MS at the intact level (top panel) and after FabRICATOR MagIC with 5 mM TCEP (bottom panel). New peaks detectable by the subunit analysis are highlighted. Data was obtained in collaboration with Symphogen A/S (Ballerup, Denmark).

Glycan analysis using FabRICATOR MagIC. Analysis of Fc glycosylation by middle-level LC-MS using FabRICATOR MagIC (top) and quantification of glycoforms from a deconvoluted mass spectrum of the Fc/2 fragment (bottom).

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| PRODUCT | DESCRIPTION | ID |
|----------------------------|-------------------------------|------------|
| FabRICATOR MagIC, 2 ml | Digests up to 24 × 200 µg IgG | A0-FRM-024 |
| FabRICATOR MagIC, 4 × 2 ml | Digests up to 96 × 200 µg IgG | A0-FRM-096 |

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
-  Human IgG1
-  KSCDKT / HTCPPC
-  O/N reaction
-  No reducing agents or co-factors needed

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Thermo Scientific™ CaptureSelect™ resin from Thermo Fisher Scientific Inc. and its subsidiaries.
Thermo Scientific and CaptureSelect are trademarks of Thermo Fisher Scientific Inc. and its subsidiaries.

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| Enzyme | FabALACTICA | GingisKHAN | Papain | Lys-C |
|---------------------|-------------------------------------|--------------------|--|--|
| Digestion site | | | | |
| Specificity | IgG-specific/ One digestion site | One digestion site | Unspecific/ Several digestion sites | Unspecific/ Several digestion sites |
| Selectivity | Human IgG1 | Human IgG1 | Several proteins | Several proteins |
| Reducing conditions | No | Yes, 2 mM cysteine | Yes | No |
| Reaction time | O/N (16-18 h) | 1 h | 1-24 h | 1-24 h |

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
Separation of intact Fab and Fc fragments by RP-HPLC after FabALACTICA digestion. The digestion was performed O/N at 37°C.

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
Paired glycan analysis of trastuzumab Fc fragments. Trastuzumab was digested using FabALACTICA O/N at 37°C, and intact Fc and Fab fragments were studied using LC-MS.

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
2D-NMR analysis following FabALACTICA digestion to evaluate the impact of oxidation on infliximab. 1H-13C methyl HSQC of the combined Fab and Fc fragments of infliximab. Originator (orange) and oxidized originator sample (green). The spectra were acquired at 310 K using an 800 MHz Bruker Avance at the Swedish NMR Center. Data was obtained in collaboration with SARomics Biostructures (Lund, Sweden).

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
Preparation of pure Fabs from three monoclonal antibodies using the FabALACTICA Fab kit. Adalimumab, trastuzumab and cetuximab digested by Immobilized FabALACTICA. Pure Fab fragments were obtained in a high yield from all three mAbs using the FabALACTICA Fab kit.

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- 


Mouse IgG2a and IgG3,
human IgG1-4, IgG of some
classes from monkey, rabbit
and sheep
- 

CPAPNLLG / GPSVF
- 

2h reaction
- 

No need for reducing
agents or co-factors

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Thermo Scientific™ CaptureSelect™ resin from Thermo Fisher Scientific Inc. and its subsidiaries.

Thermo Scientific and CaptureSelect are trademarks of Thermo Fisher Scientific Inc. and its subsidiaries.

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
Digestion of mouse IgG2a using FabRICATOR Z and FabRICATOR.

The F(ab')₂ fragment is detected at ~110kDa and the Fc fragments at ~30kDa. The enzyme is detected at ~37kDa.

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
Schematic overview of FragIT Z and FragIT Z kit. Generation and purification of F(ab')₂ and Fc fragments from mouse IgG2a and IgG3 using FragIT Z and FragIT Z kit.

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
-  Human IgG and IgG from mouse, rat, goat, sheep and rabbit.
-  KTHT / CPPCPAP (human IgG1)
-  60 min reaction
-  Requires reducing conditions

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

Generating Fab and Fc Fragments using FabULOUS. Left: SDS-PAGE analysis of human IgG1 and IgG2 and mouse IgG1 and IgG2a digested in the hinge region using FabULOUS. The generated fragments were analyzed with SDS-PAGE under reducing conditions. Right: SDS-PAGE analysis of rat IgG2b digested in the hinge region using FabULOUS at different reducing conditions. The generated fragments were buffer exchanged and analyzed with SDS-PAGE under non-reducing conditions.

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Generation and purification of Fab fragments from mouse IgG1. The top panel shows the intact monoclonal mouse IgG1 antibody, the middle panel shows the analysis of the fragments after FabULOUS digestion, and the bottom panel shows the eluted Fab fragments from the CaptureSelect™ LC-Kappa (mur) column.

FabULOUS™

| PRODUCT | DESCRIPTION | ID |
|----------------------|------------------|------------|
| FabULOUS, 2000 units | Digests 2 mg IgG | A0-PU1-020 |

FabULOUS™ Fab kit

| PRODUCT | DESCRIPTION | ID |
|-------------------------|-------------------------------------|------------|
| FabULOUS Fab kit, mouse | Digests and purifies 2 mg mouse IgG | A1-PFK-020 |

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-  Human IgG1
-  KSCDK / THTCPPCP
-  60 min reaction
-  Requires mild reducing conditions. 2mM cysteine is included

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

Specific digestion of human IgG1 using GingisKHAN. LC separation of Fab and Fc fragments from a selection of therapeutic antibodies (trastuzumab and adalimumab) and a Fc-fusion protein (etanercept) after GingisKHAN digestion.

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SDS-PAGE analysis under non-reducing conditions of the human IgG1 trastuzumab digested by GingisKHAN. After digestion, the Fab fragments were captured in the affinity spin column and could easily be eluted.

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
- 


Human IgG1-4, Fc-fusion proteins, IgG from mouse, rabbit, rat, monkey, sheep, goat, cow and horse
- 

Hydrolyzes the β 1,4-linkage between the two innermost GlcNAc residues
- 

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

LC-MS characterization of Fc glycan hydrolysis by GlycINATOR. Native and GlycINATOR-treated trastuzumab was digested with FabRICATOR to generate Fc/2 fragments before LC-MS analysis.

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
Study of ADCC by native and deglycosylated trastuzumab. The effect of trastuzumab was analyzed using iLite ADCC Effector (V) Assay Ready Cells (BM4001) and iLite ADCC Target HER2 (+) Assay Ready Cells (BM4011). Effector and Target cells were diluted and mixed with the antibody solutions. Incubation was performed at 37°C and 5% CO₂. Luminescence reading was performed after 6 h.

Immobilized GlycINATOR®

| PRODUCT | DESCRIPTION | ID |
|---|--------------------------------|-------------|
| Immobilized GlycINATOR, Microspin 2 × 0.5 mg | Deglycosylates 2 × 0.5 mg IgG | A0-GL6-010 |
| Immobilized GlycINATOR, Microspin 5 × 0.5 mg | Deglycosylates 5 × 0.5 mg IgG | A0-GL6-025 |
| Immobilized GlycINATOR, Microspin 10 × 0.5 mg | Deglycosylates 10 × 0.5 mg IgG | A0-GL6-050 |
| Immobilized GlycINATOR, Midispin 1-10 mg | Deglycosylates 1-10 mg IgG | A0-GL6-100 |
| Immobilized GlycINATOR, Maxispin 10-100 mg | Deglycosylates 10-100 mg IgG | A0-GL6-1000 |

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-  Human IgG1-4, Fc-fusion proteins, IgG from mouse, rabbit, rat, monkey, sheep, goat, cow and horse
-  Hydrolyzes the β 1,4-linkage between the two innermost GlcNAc residues
-  30 min reaction
-  Requires native IgG fold

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| | GLYCINATOR | IGGZERO | PNGASE F | ENDO H |
|-----------------------|-------------------|---|--|---|
| Digestion site | | | | |
| Works on native IgG | Yes | Yes | Yes | Yes/No |
| pH optimum | 7.4 | 7.4 | 7.5 | 5-6 |
| Reaction time | 30 min | 30 min | 15 min + | 1-16 h |
| IgG specific | Yes* | Yes | No | No |
| Glycoform specificity | All Fc glycoforms | Complex Fc glycoforms, limited activity on high-mannose and hybrid-type glycans | All glycoforms except core α 1,3 fucosylation | High-mannose and some hybrid-type glycans |

* GlycINATOR has one more known substrate, α -1-acid glycoprotein.

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
References

Henninot, A. et al. *Characterization of monoclonal antibodies by a fast and easy LC- MS ToF analysis on culture supernatant.* Anal Biochem 1-8 (2015).

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
- 

Human IgG1-4, IgG from mouse, rabbit, rat, monkey, sheep, goat, cow and horse
- 

Conjugation occurs at the Fc N-glycan sites
- 

2-3 day protocol
- 

Alexa Fluor®, biotin, DFO, MMAE, PNU, Azide activation

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
Analysis of Fc glycan conjugation using GlyCLICK. RP-HPLC analysis of unmodified panitumumab and panitumumab conjugated with desferrioxamine (DFO), Alexa Fluor® 488, biotin or monomethyl auristatin E (MMAE) after digestion with FabRICATOR.


1. SiteClick™ is provided under an intellectual property license from Life Technologies Corporation. The trademark SiteClick™ is the property of Life Technologies Corporation.
2. Linker-toxin Payloads included in GlyCLICK® are provided under an intellectual property license from Glykos Finland Ltd.

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Affinity analysis of native and conjugated trastuzumab. Anti-human IgG (Fc) was used as the capturing molecule for trastuzumab. Native (control), Alexa Fluor® 488 (GlyCLICK, DOL=2), DyLight® 488 (random, DOL=10). HER2 was injected in a range to ensure sufficient curvature. All data were fitted against a 1:1 mathematical model.

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

Tumor uptake

Biodistribution ex vivo

Analysis of PET/CT imaging of SK-OV-3 tumor-bearing mice. Mean tumor uptake (left panel) over a time interval of 0-168 hours post injection, and ex vivo biodistribution of ^{89}Zr -DFO-trastuzumab in major organs (right panel). Data was obtained in collaboration with Minerva Imaging (Copenhagen, Denmark).

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


- 
Human IgG1, IgG2 and IgG4
- 
Transglycosylation occurs at the Fc N-glycan sites
- 
3 hour workflow
- 
Glycoform included, optional defucosylation using the Afucosylated kit

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TransGLYCIT™

| PRODUCT | DESCRIPTION | ID |
|------------------------|--|------------|
| TransGLYCIT G0, 1 mg | Generates 1 mg human IgG with the G0 glycoform | T1-G0F-010 |
| TransGLYCIT G1, 1 mg | Generates 1 mg human IgG with the G1 glycoform | T1-G1F-010 |
| TransGLYCIT G2, 1 mg | Generates 1 mg human IgG with the G2 glycoform | T1-G2F-010 |
| TransGLYCIT G2S2, 1 mg | Generates 1 mg human IgG with the G2S2 glycoform | T1-S2F-010 |

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
Transglycosylation using TransGLYCIT. Deconvoluted mass spectra of the Fc/2 fragment of native trastuzumab (top) and after transglycosylation with TransGLYCIT (left panel) or TransGLYCIT Afucosylated (right panel). The mAb was digested with FABRICATOR and the subunits were analyzed by reversed-phase LC-MS on a Waters™ BioAccord™ system equipped with a Waters™ BioResolve™ RP mAb column (2.1 × 50 mm).

TransGLYCIT™ Afucosylated

| PRODUCT | DESCRIPTION | ID |
|-------------------------------------|---|------------|
| TransGLYCIT G0 Afucosylated, 1 mg | Generates 1 mg afucosylated human IgG with the G0 glycoform | T1-G0A-010 |
| TransGLYCIT G1 Afucosylated, 1 mg | Generates 1 mg afucosylated human IgG with the G1 glycoform | T1-G1A-010 |
| TransGLYCIT G2 Afucosylated, 1 mg | Generates 1 mg afucosylated human IgG with the G2 glycoform | T1-G2A-010 |
| TransGLYCIT G2S2 Afucosylated, 1 mg | Generates 1 mg afucosylated human IgG with the G2S2 glycoform | T1-S2A-010 |

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-  Digests any peptide or protein containing Arg. Specific for Arg-X motifs
-  Digests C-terminally of Arg residues
-  60 min reaction
-  Active in 6M urea and 0.1% SDS

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
Digestion of oxidized β-chain of insulin with GingisREX and Arg-C.

The digestion was performed O/N at 37°C, with an enzyme to substrate ratio of 1:20 (w/w), 20mM cysteine in buffers at pH 7.4 (GingisREX) and pH 7.6 (Arg-C). The peptides were separated on RP-HPLC.


| PEPTIDE NO. | AMINO ACID SEQUENCE |
|----------------|-------------------------------|
| Intact protein | FVNQHLGSHLVEALYLVCGERGFFYTPKA |
| 1 | GFFYTPKA |
| 2 | FVNQHLGSHLVEALYLVCGER |
| 3 | GFFYTPK |
| 4 | FVNQHLGSH |
| 5 | LVEALYLVCGER + Na |

Sequences of oxidized insulin β-chain digested by GingisREX or Arg-C. Green indicates arginine residues and red indicates lysine residues.

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-  N-linked glycans on glycoproteins
-  Hydrolyzes the glycosidic bond between N-glycans and Asn
-  15+ min denaturing or 1+h native reaction
-  RapiGest™ included in the Denaturing kit

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
Immobilized PNGase F

| PRODUCT | DESCRIPTION | ID |
|--|---------------------------------------|------------|
| Immobilized PNGase F, Microspin 5×0.2mg | Deglycosylates 5×0.2mg glycoproteins | G1-PF6-010 |
| Immobilized PNGase F, Microspin 10×0.2mg | Deglycosylates 10×0.2mg glycoproteins | G1-PF6-020 |

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
Removal of N-glycans under native conditions. TIC chromatogram (left) and deconvoluted mass spectra (right) of the Fc/2 fragment of intact trastuzumab (top panel) and trastuzumab treated with Immobilized PNGase F (bottom panel) or PNGase F in solution (middle panel). The antibody samples were analyzed after digestion with FabRICATOR.

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Key characteristics and deconvoluted mass spectra or SDS-PAGE assay of a selection of glycoproteins treated with Immobilized PNGase F. Abatacept was analyzed by reversed-phase LC-MS and RNase B and cetuximab were analyzed on SDS-PAGE. The mass shifts in the deconvoluted mass spectra or on the gels show the successful removal of N-glycans from the various glycoprotein substrates.

Immobilized PNGase F Denaturing


| PRODUCT | DESCRIPTION | ID |
|---|---|------------|
| Immobilized PNGase F Denaturing, Microspin 5×0.2mg | Deglycosylates 5×0.2mg glycoprotein. 5×1 mg <i>RapiGest</i> [™] SF is included | G2-PDK-010 |
| Immobilized PNGase F Denaturing, Microspin 10×0.2mg | Deglycosylates 10×0.2mg glycoprotein. 10×1 mg <i>RapiGest</i> [™] SF is included | G2-PDK-020 |

RapiGest[™] SF Surfactant is included in Immobilized PNGase F Denaturing. *RapiGest*[™] SF Surfactant is a trademark of Waters Corporation.

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- 

Native proteins with mucin-type O-glycosylation
- 

Digests N-terminally of O-glycosylated Ser and Thr residues
- 

2h to O/N (16-18h) reaction
- 

Desialylation using SialEXO (included) increases performance

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OpeRATOR®

| PRODUCT | DESCRIPTION | ID |
|----------------------|-----------------------------|------------|
| OpeRATOR, 2000 units | Digests 2 mg O-glycoprotein | G2-OP1-020 |

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
Specific digestion N-terminally of the O-glycosylation site. The reduced fragments were separated on a reversed-phase C4 column followed by ESI-QTOF MS detection. The EPO protein carrying one core 1 O-glycan was hydrolyzed at a single specific site N-terminally of the O-glycosylated serine using OpeRATOR.

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

In-depth analysis of the O-glycosylation sites of etanercept originator. HILIC-MS analysis (top panel) of etanercept peptides generated by digestion with OpeRATOR. The middle panel shows the peptide map of the O-glycosylation sites of etanercept based on the results from the OpeRATOR-based workflow. The identified O-glycosylation sites are marked in yellow. Two sites (marked with an asterisk) could only be inferred from the OpeRATOR digestion pattern without direct identification of a peptide containing the glycosylated amino acid. The bottom panel shows a RP-MS analysis of etanercept digested with trypsin. Data was obtained in collaboration with Thermo Fisher Scientific (UK).

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-  O-glycans on glycoproteins
-  Hydrolyzes Core 1 disaccharides
-  2-4 h reaction
-  Requires removal of terminal sialic acids (SialEXO included)

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
Deglycosylation of the TNFR subunit of etanercept from three commercially available preparations.

Deconvoluted mass spectra. Left panel: The TNFR subunit carrying both N- and O-glycans is too complex for all variants to be resolved. Middle panel: Treatment with PNGase F for removal of N-glycans enabled assessment of the O-glycan composition. C1S1 = core 1 monosialylated, C1S2 = core 1 disialylated. Right panel: Complete glycan removal (using OglyZOR, SialEXO, PNGaseF) provided information of the primary structure of the TNFR part. N-terminal truncation, -L (2-256) and -LP (3-256) are visualized. A fraction of the deglycosylated subunits revealed mass shifts of +203 Da (HexNAc), indicating presence of incomplete core 1, not removed by OglyZOR.

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
Comparison of the enzymatic activities of OglyZOR and SialEXO to commercially available endoglycosidases and sialidases.

All incubations (4 h) were performed according to the manufacturer's instructions, and the samples were analyzed by SDS-PAGE.

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- Sialic acids on N- and O-linked glycans
- α 2-3, α 2-6 and α 2-8-linked sialic acids
- 30 min – 2 h reaction
- Requires no co-factors

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
Analysis of released sialic acids. HILIC analysis of released glycans with α 2-3-linked (left) or α 2-6-linked (right) sialic acids after incubation with SialEXO 23 or SialEXO respectively. Sialylated glycan structures are shaded in purple.

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
Unmasking charge variants. The analytical workflow improved the separation and enabled charge variant analysis of glycoproteins. Etanercept is difficult to analyze without first removing the sialic acids, but using Immobilized SialEXO, separated peaks were obtained. Data was obtained in collaboration with ProteinSimple (San Jose, CA, USA).

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-  Fucose on N- and O-glycosylated proteins
-  Hydrolyzes α1-2,3,4-linked fucose
-  1-2h reaction
-  No need for co-factors

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FucosEXO™


| PRODUCT | DESCRIPTION | ID |
|----------------------|---------------------------------|------------|
| FucosEXO, 2000 units | Defucosylates 2 mg glycoprotein | G1-FM1-020 |

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FucosEXO substrate specificity. The substrate specificity of FucosEXO was analyzed on equal molar amounts of synthetic oligosaccharides; α 1-2 (2'-fucosyllactose), α 1-3 (3-fucosyllactose), α 1-4 (Lewis a) and α 1-6 (α 1-6 fucosylated chitobiose). Substrates were incubated with FucosEXO for 30 min at 37°C and the amount of released fucose was measured spectrophotometrically using an L-Fucose Assay Kit (Megazyme).

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
Defucosylation using Immobilized FucosEXO. Glycoengineered TNFR with core 1 and core 2 O-glycans decorated with both α 1-2 and α 1-3 fucose was incubated with bovine kidney fucosidase (1 h, 37°C) or Immobilized FucosEXO (1 h, RT). The processed protein was analyzed by reversed-phase LC-MS on a Waters™ BioAccord™ system equipped with a Waters™ BioResolve™ RP mAb column (2.1 × 50 mm).

Glycoengineered TNFR (PNGase F treated)

Immobilized FucosEXO™

| PRODUCT | DESCRIPTION | ID |
|---|--|------------|
| Immobilized FucosEXO, Microspin 5 × 0.5 mg | Defucosylates 5 × 0.5 mg glycoprotein | G1-FM6-025 |
| Immobilized FucosEXO, Microspin 10 × 0.5 mg | Defucosylates 10 × 0.5 mg glycoprotein | G1-FM6-050 |

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-  Galactose on N- and O-glycosylated proteins
-  Hydrolyzes β1-3 and β1-4-linked galactose
-  2h reaction
-  No co-factors needed

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


Released glycan analysis. HILIC-FLD UHPLC chromatograms of a 2-AB labeled glycan library analyzed undigested (top), after treatment with SialEXO (middle) or SialEXO and GalactEXO (bottom). Analyzed on Thermo Scientific Vanquish Duo UHPLC system equipped with Thermo Scientific Accurcore 150 Amide HILIC column (2.1 × 150 mm).

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Degalactosylation of trastuzumab. Deconvoluted mass spectra of the Fc/2 fragment of trastuzumab treated with immobilized GalactEXO (bottom panel). The mAb was digested with FabRICATOR, separated by RP-HPLC (Waters™ BioResolve™ RP mAb (2.1 × 50 mm)) and analyzed by ESI-Q-TOF MS (Bruker Impact II).

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-  α-linked GalNAc residues on O-glycoproteins
-  Tn antigen and α1-3-linked terminal GalNAcs
-  4 h reaction
-  No need for co-factors

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
GalNAcEXO activity on etanercept with truncated O-glycans. Deconvoluted mass spectra of partially deglycosylated etanercept (top) and after treatment with GalNAcEXO (O/N at 37°C, bottom). Samples were FabRICATOR-digested, separated by RP-HPLC (Waters™ BioResolve™ RP mAb, 2.7 μm, 2.1 × 100 mm) and analyzed by ESI-Q-TOF MS (Bruker Impact II).

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
Total deglycosylation of the recombinant C1 inhibitor. The C1 inhibitor was analyzed by LC (Waters™ ACQUITY UPLC Protein BEH C4 column, 1.7 μm, 2.1 × 100 mm) and ESI-Q-TOF mass spectrometry (Bruker Impact II). In its intact form (left), the protein is highly complex, and deconvolution of the mass spectra only yielded noise. After removal of N- and core 1 O-glycans, the Tn antigens remained (middle) and could efficiently be removed by incubation with GalNAcEXO (O/N at 37°C, right).

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-  Mucin-type O-glycosylated proteins and peptides
-  Binds to glycoproteins and peptides carrying mucin-type O-glycans.
-  2h workflow
-  Desialylation using SialEXO (included) increases the binding

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GlycOCATCH®

| PRODUCT | DESCRIPTION | ID |
|------------|-------------------------------------|------------|
| GlycOCATCH | Enrichment of 0.2 mg O-glycoprotein | G3-OC6-002 |

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
Affinity purification using GlycOCATCH.

Selective purification of an O-glycosylated peptide (glycodrosocin) carrying a core 1 O-glycan.

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
Enrichment of O-glycoproteins from human serum. Human serum was treated with SialEXO and incubated with the GlycOCATCH resin. The eluted proteins were denatured, reduced, alkylated and trypsin digested. Peptides were analyzed using RP-LC-MS/MS, and searched against the Swiss Prot database. Identified O-glycosylated proteins and other proteins with > 12 matching peptides and a MASCOT score > 200 are listed.

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- 
Detects the FabRICATOR (IdeS) enzyme
- 
Source: goat
- 
Polyclonal
- 
Concentrations ranging from 0.5 mg/ml to 4 mg/ml

**Anti-FabRICATOR®
Affinity Purified**

**Anti-FabRICATOR®
Affinity Purified Biotin Conjugated**


**Anti-FabRICATOR®
Protein G Purified**

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Performance of Anti-FabRICATOR Affinity Purified. To display the difference in performance between the Anti-FabRICATOR Protein G Purified and the Anti-FabRICATOR Affinity Purified, wells were coated with either Anti-FabRICATOR Protein G Purified or Anti-FabRICATOR Affinity Purified at 2 µg/ml. The FabRICATOR enzyme was added in dilutions from 1 µg/ml to 1 ng/ml. Anti-FabRICATOR Affinity Purified Biotin Conjugated was then added at 1 µg/ml. HRP-Streptavidin was added in a 1:2000 dilution, followed by visualization using ABTS.

Anti-FabRICATOR®

| PRODUCT | DESCRIPTION | ID |
|--|--|------------|
| Anti-FabRICATOR Affinity Purified | Frozen solution of 0.1 ml, concentration 1 mg/ml | A3-AF2-010 |
| Anti-FabRICATOR Affinity Purified, Biotin Conjugated | Frozen solution of 0.1 ml, concentration 0.5 mg/ml | A3-AF3-005 |
| Anti-FabRICATOR Protein G Purified | Frozen solution of 0.1 ml, concentration 4 mg/ml | A3-AF1-010 |

Legal and Disclaimers

All rights reserved. Genovis products may be covered by one or more patents, trademarks and copyrights owned or controlled by Genovis AB. For more information about commercial rights, please contact the Genovis team at licensing@genovis.com.

Genovis products are intended for research use only. They are not intended to be used for therapeutic or diagnostic purposes in humans or animals. All goods and services are sold subject to Genovis' General Terms and Conditions of Sale.

FabRICATOR®

This product is provided under an exclusive world-wide intellectual property license from Hansa Biopharma AB derived from international publication WO03051914, including granted US Patent No US 7,666,582 and granted European Patent No 1458861. The license encompasses IdeS from *Streptococcus pyogenes* for biotechnical industrial applications which are neither therapeutic nor diagnostic, other than the following exception which is included within the license: digesting IgG *in vitro* in clinical samples for diagnostic purposes.

POROS® included in FabRICATOR®HPLC

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research or for analytical testing during manufacturing, all conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components

for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; and/or (c) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Thermo Fisher Scientific, 5826 Newton Drive, Carlsbad, CA 92008 USA or outlicensing@thermofisher.com. The trademark POROS® trademark is the property of Life Technologies Corporation.

FabRICATOR® Magic

The trademarks Thermo Scientific™, KingFisher™ and BindIt™ are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

SiteClick™ included in GlyCLICK®

This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than

as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@thermofisher.com. The trademark SiteClick™ is the property of Life Technologies Corporation.

Linker-toxin Payloads included in GlyCLICK®

This product is provided under an intellectual property license from Glykos Finland Ltd. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer. The buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; or (b) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research.

CaptureSelect™

Thermo Scientific™ CaptureSelect™ resin from Thermo Fisher Scientific Inc. and its subsidiaries. Thermo Scientific and CaptureSelect are trademarks of Thermo Fisher Scientific Inc. and its subsidiaries.

TransGLYCIT™

TransGLYCIT contains a technology licensed from GlycoT Therapeutics, LLC.

Immobilized PNGase F Denaturing

RapiGest™ SF Surfactant from Waters Corporation is included in Immobilized PNGase F Denaturing. *RapiGest™* is a trademark of Waters Corporation.

© Genovis AB

The SmartEnzymes portfolio consists of enzymatic products and technologies designed to improve the efficacy and throughput in analytical or preparative workflows for complex biopharmaceuticals such as antibodies, Fc-fusion proteins, complex glycoproteins and antibody conjugates.

FabRICATOR®
Below hinge
digestion of IgG

FabALACTICA®
Above hinge
digestion of human IgG1

FabRICATOR® Z
Below hinge
digestion of mouse IgG

FabULOUS™
Above hinge
digestion of IgG

GingisKHAN®
Above hinge
digestion of human IgG1

GlycINATOR®
Hydrolysis of
IgG Fc glycans

IgGZERO®
Hydrolysis of
complex Fc glycans

GlyCLICK®
Site-specific
conjugation of native IgG

TransGLYCIT™
Transglycosylation
of IgG

GingisREX®
Arginine-specific
protein digestion

**Immobilized
PNGase F**
Hydrolysis of N-glycans

OpeRATOR®
O-glycan-specific
protein digestion

OglyZOR®
Hydrolysis of
core-1 O-glycans

SialEXO®
Hydrolysis of
sialic acids

FucoseEXO™
Hydrolysis of
 α 1-2,3,4 fucose

GalactEXO™
Hydrolysis of β 1-3,4-
linked galactose

GalNAcEXO™
Hydrolysis of
 α -linked GalNAcs

GlycOCATCH®
Enrichment of
O-glycoproteins

Anti-FabRICATOR®
Detection of the
FabRICATOR enzyme

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Shigematsu & Co., Ltd.