





光分解性修飾 / Photocleavable Modifiers

製品番号	製品名 / 化学名	構造式	製品情報
2130	PC 5'-Amino-Modifier-CE Phosphoramidite 1-[2-Nitro-5-(6-trifluoroacetylcaproamidomethyl)phenyl]-ethyl-[2-cyano-ethyl-(N,N-diisopropyl)]-phosphoramidite		【パッキング】 100μmol 250mg 【劇物】
2122	PC 5'-Biotin-CE Phosphoramidite 1-[2-Nitro-5-(6-(N-(4,4'-dimethoxytrityl))-biotinamidocaproamidomethyl)phenyl]-ethyl-[2-cyanoethyl-(N,N-diisopropyl)]-phosphoramidite		【パッキング】 50μmol 100μmol 250mg 【劇物】
2131	PC Spacer-CE Phosphoramidite 1-[2-Nitro-5-(6-(4,4'-dimethoxytrityloxy)butyramidomethyl)phenyl]-ethyl-[2-cyanoethyl-(N,N-diisopropyl)]-phosphoramidite		【パッキング】 100μmol 250mg 【劇物】
2066	PC Linker-CE Phosphoramidite 3-(4,4'-Dimethoxytrityl)-1-(2-nitrophenyl)-propane-1,3-diol-[2-cyanoethyl-(N,N-diisopropyl)]-phosphoramidite		【パッキング】 100μmol 250mg 【劇物】